

General Hospital

Correlation between HER2 and ALK status in circulating tumor Cells (CTCs) and tissue of breast and lung cancer patients

Chassidy Johnson¹, Ali Asgar Bhagat¹, Man Chun Leong¹, Tse Hui Lim², Yong Wei Chua², Garima Singh¹, Andrew Wu¹, Yifang Lee¹, Alvin Soon-Tiong Lim², Tony Kiat Hon Lim², Darren Wan-Teck Lim³, Yoon-Sim Yap³



National Cancer Centre Singapore SingHealth

¹Clearbridge Biomedics Pte Ltd, Singapore; ²Department of Pathology, Singapore General Hospital, Singapore; ³National Cancer Centre Singapore, Singapore

Introduction

Evaluation of CTCs offers a minimally invasive means to access tumor heterogeneity and disease evolution in response to selective pressures of targeted therapies in real-time, making them useful tools to guide treatment decisions. Here we present data on the concordance of HER2 and ALK status between tissues and CTCs, and the utility of these biomarkers on CTCs to examine heterogeneity of disease and potentially prognosticate outcome in breast and NSCLC patients.

Methods

CTC Enrichment by ClearCell® FX and FISH analysis



Figure 1. Workflow: Whole breast blood (n=25) and cancer NSCLC (n=27) patients was collected in EDTA tubes. 7.5ml of blood was RBC lysed, and the nucleated cells were CTC for processed enrichment using the ClearCell FX system. O† enriched were CTCs fixed, cytospun onto slides and glass for FISH processed analysis.

Fluorescence in-situ hybridization (FISH) on CTCs: Samples were fixed with ethanol series, denatured, and hybridized with FISH probes specific for HER2 amplification and ALK rearrangement. FISH scoring was performed under an epifluorescence microscope at 1000× magnification by certified cytogeneticists blinded to the previously defined HER2 and ALK profiles. **ClearCell[®] FX system enrichment principle:** Patient blood (sample) and sheath fluid enter the microfluidic chip, creating a density gradient (X). Curvature in the tube creates shear and lift forces causing cell size-dependent migration of cells across the established gradient density (Y). The larger cells (CTCs) are optimally separated from smaller blood cells at a distance of 1 dean cycle where the outlets collecting enriched CTCs (sample) and blood cells (waste) are located (Z).

Table 1. Patient Demo cohort	ographics of NSCLC	Table 2. Patient Demographics of Breast Can cohort			
Patient characteristics	Cases (%), N= 27	Patient characteristics	Cases (%), N= 25		
Age, years	32–76	Age, years	36-66		
Gender		Gender			
Male	16 (59.3%)	Male	0 (0%)		
Female	11 (40.7%)	Female	25 (100%)		
Smoking history		Clinical staging			
Non-smoker	16 (59.3%)	111	8 (32%)		
Smoker	5 (18.5%)	IV	17 (68%)		
Ex-smoker	5 (18.5%)	Histological subtype			
No info	1 (3.7%)	HER2-positive	16 (64%)		
Clinical staging		HER2-negative	9 (36%)		
IB	1 (3.7%)				
IIIA	1 (3.7%)				
IIIB	2 (7.4%)				
IV	23 (85.2%)				
Histological subtype					
ALK-positive	14 (51.9%)				
ALK-negative	12 (44.4%)				

Therapy Naïve for *ALK*--targeted TKI 27/27 (100%)

Results

cer





Strong correlation between *HER2* and *ALK* status between tissue and CTCs

4. [• ALK		Tumor tissue		В.			Tumor tissue	
			ALK (+) †	ALK (-)				HER2 (+) [†]	HER2 (-)
	СТС	ALK (+)*	12	0		CTC	HER2 (+) *	10	1
	010	ALK (-)	0	12			HER2 (-)	6	8

62.5% Sensitivity/ 88.9% Specificity (Based on ROC curve) 100% Sensitivity/Specificity (Based Off ROC Curve) *+ ALK*+ tissue cut point >15% cell positivity *+ HER2*+ tissue cut point IHC 3+ and IHC2+ FISH confirmation * ALK+ CTC based on empirical cut-off >2 CTCs/~2 ml blood * HER2+ CTC empirical cut-off of > 3 cells/7.5 ml blood

Figure 3. (A) Patients were defined as *ALK*(+) using >15% of tumor cells in FFPE tissue and >2 CTCs positive for ALK rearrangement in tissue and CTCs, respectively. ALK rearrangements in CTCs were identified in 12/12 (100%) ALK(+) patients and 0/12 (0%) of ALK(-) patients. Based on ROC curve analysis, sensitivity = 100% and specificity = 100%. (B) Cells identified with HER2 gene amplification and CEP17 polysomy are scored as "positive cells". HER2(+) CTCs were identified in 10/16 (62.5%) HER2(+) patients. HER2 gene-amplified CTCs are identified in 7/16 (43.8%) HER2(+) patients, range (0-12 cells/7.5ml blood, N=16 patients). 0/9 (0%) HER2(-) patients were identified with HER2 geneamplified CTCs, although CEP17 polysomy cells are observed in 5/9 (55.6%) patients. Based on ROC curve analysis, 62.5% sensitivity and 88.9% specificity.

Association between number of *HER2*-positive CTCs at baseline and trastuzumab response



Figure 4. Kaplan-Meier survival curve of HER2-positive tumor patients undergoing trastuzumab treatment. High *HER2*-positive cell cut-off defined as \geq 5 cells from 7.5 ml analyzed blood. Median PFS for *HER2* high (> 5 cells) vs HER2 low (<5 cells) is 10.5 vs 4.5 months; HR 1.934.



Figure 5. Representative case studies where HER2 status in CTCs was evaluated at baseline, on trastuzumab treatment (on-TRX) and at progression for patients 12 (left) and 13 (right) are shown. CTCs defined as HER2 amplification (gold) and polysomy 17 (grey) are represented as percent of total CTCs detected and actual cell numbers are indicated. The trend in ratio of *HER2* amplified CTCs is summarized. * Sample at progression for Patient 12 was not available.

The frequency and ratio of CTCs with HER2 amplification changes during treatment, highlighting the evolution of disease under selective pressures of targeted therapy. CTCs may provide a useful tool to evaluate these changes, capture heterogeneity of disease during treatment and potentially help guide treatment decisions.

Baseline (0 mo			
6			
1F1R1G			
1F1R1G			

Clinical response to Crizotinib correlates with ALK positive CTCs

Figure 6. A case study from a never smoker, male, diagnosed with NSCLC and no accessible tissue for ALK FISH testing is shown. CT Scans and blood draws to examine ALK rearrangement patterns in CTCs was performed at baseline (0 months), partial response (3 months) and progression (5 months). Shown are the CT scans illustrating presence of metastatic tumor(s) on the liver (yellow circles), number of ALK+ CTCs, rearrangement patterns detected and representative FISH images for the respective time points. The rearrangement pattern in the image is indicated.

The results presented contribute to the growing body of evidence regarding the utility of using CTCs as tools to access the heterogeneity and evolution of cancer when under selective pressures of targeted therapy. The strong concordance of HER2 and ALK status between tissue and CTCs in breast and NSCLC warrants further investigation into using CTCs enriched using the ClearCell[®] FX system to potentially help guide treatment decisions.



HER2 status in CTCs: Analysis of heterogeneity and disease evolution in patients on trastuzumab treatment

Partial Response (3 months) Progression (5 months) 1F1R1G 1F1R1G, 2R2G, 1F1R 2R2G **1F1R**

Conclusions