# Whole transcriptomics analyses of mimicking Circulating Tumor Cells (CTCs) by single-cell RNA sequencing (scRNAseq)

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### Purposes

Ultra-precision medicine is emerging as a cornerstone of cancer biology, as it aims a providing a better comprehension of cancer cells and their clinical response to immunotherapy. Consistently, single-cell sequencing of CTCs is a powerful tool to further decipher tumor plasticity and to identify new pharmacological targets influencing clinical outcome and response to treatment. The aim of this study is to validate an original workflow allowing the isolation at the single-cell level of mimicking CTCs without interfering with singlecell RNAseq analysis.

### Methodology

eliminated with a lysis buffer. The residual white blood cells and mCTC's were loaded in the ClearCell CTChip FR1. mCTC The technological advances in microfluidic systems and isolation technologies have resulted enrichment is based on the Dean Flow Fractionation principle. The smaller hematological cells [8 – 15 µm] are affected by the in the enriched extraction of mimicking CTCs from healthy whole blood samples. In the Dean Drag and migrate to waste outlet whereas the mCTC's [>15 µm] migrate to the upper outlet. present study, the ClearCell Fx (Biolidics Limited) was used as a label-free microfluidic system for enrichment of wholly intact CTCs, while the cellenONE F1.4 system (Cellenion) Experiment design Condition #1 Physical integrity of the cell mCTC enriched with the Clearcell was used to isolate single CTCs. The latter allowed high-throughput automated isolation and (Only mCTC) Fluorescent FX from a total blood sample shape after sorting mCTC dispensing of single CTCs in 96-well plates containing cell lysis buffer. From these 96-well Conservation of cell circularity 23 20 3 (86.9) plates, scRNA libraries were prepared with the NebNext Single Cell/Low Input kit (New #1 \_\_\_\_\_\_#2 and elongation after th England Biolabs) and sequenced with paired-end sequencing (2x75 bp) on the NextSeq ClearCell and CellenOne steps ClearCell Fx Mean Area: 367 CellenOne 500/550 from Illumina. FASTQ files, obtained by the demultiplexing of Base-Call Files (BCL) Mean circularity: 1,04 Mean elongation: 1.62 were used for reads alignment and gene annotation based on the Hg19 (GRCh37.p5) CellenOne Achievement of a recovery rate Biolidics 🕺 cellenio genome reference. The mapped and annotated reads were stored in BAM files, processed of 86.9% of isolated mCTCs mixed with WBC using the Gfold (Generalized FoldChange) algorithm to produced normalized read counts. Those were then uploaded in a web-based platform named ASAP (Automated Single-cell Analysis Pipeline) to generate PCA analyses in order to distinguish cell populations and to **Figure 5:** Experimental design to evaluate the CellenOne sorting impact on mCTC morphology and recovery rate produce HEATmaps to observe the differential expression of genes between cells.

heterogeneity.

properties during embryogenesis and cancer progression.

Figure 3: Overview of the CTC isolation workflow workflow at the single cell level and comparison with a bulk of 200 mCTCs Authors would like to thank all patients, families and prescribing physicians.

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point 2 named C2 time point 3 named B. (B) Principle dysregulated transcriptomics signatures in the Bulk condition (B) compared to trypsinized cells (C1)

Figure 4 : Enrichment of Circulating Tumor Cells (CTS's) by a Label-Free Inertial microfluidic method: the red blood cells were

## Results

pathogen.

Figure 10: Log2 fold-changes of markers, expressed in breast cancer cells

Figure 12: Log2 fold-changes of transcription factors in bulk of 200 HMECs with and without the similar frozen step between Clearcell FX1 enrichment and CellenOne isolation. ERCC positive control is a pool of 50 cell lines and m\_HMEC is the transcriptomics of the mesenchymal sub-population of HMECs.

In conclusion, we propose this workflow as a standard protocol to analyze CTCs owing to its specificity and reliability at single-cell level. Such knowledge of single-cell biology may lead to the development of specific therapies to limit tumor progression and seeding of tumoral cells into secondary healthy organs by blocking newly identified targets

ZEB1, ZEB2 and twist 1 are over expressed in mesenchymal HMEC cells

### **Conclusion & Perspectives**